



POROLUX™ Revo series




The POROLUX™ Revo series is the revolution in porometry. Setting the bar in pressure step/stability method with our patent pending MP² (Multistage Pressure Process) technology, the POROLUX™ Revo series delivers the most accurate and reproducible pore size measurements, in the highest resolution.

Thanks to its cutting-edge PoreSmart software, analyzing filter media is very straightforward. Built-in intelligence and unique features such as the re-evaluation function make the instrument a pleasure to work with.

The POROLUX™ Revo's are the right instruments to characterize a wide range of materials with complex pore structures, such as polymeric membranes (flat sheets and hollow fibers), ceramic membranes, porous metals and nonwovens.

Key features

- Delivers the most accurate and reproducible pore size measurements
- Very fast and reproducible determination of the bubble point (largest pore according to ASTM F-316-03), mean flow pore size, smallest pore, pore size distribution, cumulative flow distribution & gas permeability
- 3 ways to calculate the bubble point:
 - BP dPL is the bubble point measured as a deviation from the linearity of a user-defined pressure increase ⁽²⁾
 - BP x-ml is the bubble point measured at a user-defined flow rate
 - BP pCF is the bubble point measured as a user-defined percentage of the cumulative flow
- Pressure step/stability method with our patent pending MP² technology
- With our enhanced mathematical model, the following additional results can be calculated: total pore number, total pore area and open porosity
- Automatic switch for pressure and flow sensors
- PoreSmart, our very intuitive and easy to use software, with unique features such as:
 - Re-evaluation button which allows to do a re-calculation on already performed tests
 - Built-in intelligence, e.g. algorithms for calculating bubble point, checking dry curve, wet/dry convergence, input pressure
 - Two-stage curve fitting (wet & dry curve and pore size distribution)
- Ethernet and USB connection, making remote access for installation, support and diagnosis via the internet possible.



Measuring curves and resulting parameters in Capillary Flow Porometry
BP = largest pore

POROLUX™ Revo series technical specifications

	Revo X	Revo
Technique	Gas-liquid porometry	Gas-liquid porometry
Measurement method	Pressure step/stability with patent pending MP ² technology	Pressure step/stability with patent pending MP ² technology
Max pressure	8 bar/116 psi	35 bar/500 psi
Min pore ⁽¹⁾	80 nm	13 nm
Max pore ⁽¹⁾	500 µm	500 µm
Bubble point	BP x-ml, BP pCF ^(*)	BP dPL, BP x-ml, BP pCF ^(*)
Max flow	200l/min	200l/min
Dimensions (DxWxH)	530x530x755 mm	530x530x755 mm
Weight	66 kg	70 kg

(1) depending on the wetting liquid

(2) only available on POROLUX™ Revo

(*) BP dPL is the bubble point measured as a deviation from the linearity of a user-defined pressure increase. BP x-ml is the bubble point measured at a user-defined flow rate. BP pCF is the bubble point measured as a user-defined percentage of the cumulative flow.

About us

For more than 15 years, we have been manufacturing porometers of the highest quality and reliability. All of our porometers have been completely designed and built in-house, using the latest techniques and materials. This enabled us to create an installed base of hundreds of instruments around the world, both in R&D and quality control labs in corporations, universities and research centers.

Experts in porometry

Purchasing a porometer does not only mean acquiring a device, it also guarantees lifelong support and advice from our team of porometry experts. Measuring daily in our labs allows us to gain a lot of knowledge and insights in porometry. Knowledge we like to share with our customers by, amongst others, helping to find the right settings for the samples.


Wide range of solutions

As our porometer range covers all common capillary flow techniques – including gas liquid and liquid liquid porometers – and measurement methods (pressure scan and pressure step stability), we can offer the best porometer solution for any given application.

What's more, our porometer portfolio covers the widest possible pore size range: from 2 nm up to 500 μm . Potential clients that aren't sure which porometer is right for them, are welcome to have some samples analyzed free of charge by our specialized application labs.

PoreSmart software

PoreSmart is synonym for smart porometry software, instrumental for our customers to generate key insights into their materials. With features such as easy parameter setting, built-in intelligence and presentable outputs eg. in Excel, PoreSmart is the reference in porometry software.



Global presence

Porometer has a global presence, with manufacturing sites, labs and offices in Europe, USA and China. Furthermore, we have an extensive network of highly trained distributors around the globe. This enables us to provide a first-class service anywhere in the world and to quickly follow up on any technical questions.



Part of Aptco Group

Our instruments are designed and manufactured in-house by our engineering and production department in Germany. Porometer is a brand of Aptco Technologies, a manufacturer of measurement instruments and testing equipment for academic and industrial quality control and research labs.

Aptco Technologies is part of Aptco Group, an international technology group of companies active in the distribution, manufacturing, servicing and calibration of scientific instruments and equipment for industrial, medical and academic laboratories.



Some of our customers



Porometer is a brand of
Aptco Technologies NV

Aptco Technologies NV
Nieuwe Steenweg 20A
9810 Nazareth
Belgium

T: +32 9 252 25 35
E: info@porometer.com
W: www.porometer.com