LABODEST®

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FISCHER®
technology
... a brand of i-Fischer Engineering GmbH

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Construction of the equilibrium curves

The phase equilibrium of a binary mixture is graphically shown by drawing the liquid composition of the boiling mixture (x) and the vapour composition (y) with the concentration each indicated in mol-% of the light-boiling component. Such a diagram can be made both for constant pressure and constant temperature. For technical purposes, it is preferred to use the diagram for constant pressure whereas thermodynamic calculations in general are based on a constant temperature.

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
Fig. 3
Equilibrium curves for different mixtures

- 1 Ethanol – water
- 2 Chloroform – benzene
- 3 Water – n-butanol
- 4 Nitric acid - water

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Evaporation is effected by means of an electrical immersion heater (1) made of quartz glass which is arranged concentrically in a flow heater. An even circulation of both phases with simultaneous magnetic stirring of the reflowing circulation streams in a special mixing chamber (2) ensures a quick equilibrium adjustment. Before entering the separation chamber, the vapour-liquid mixture passes a lengthened contact path (3) which guarantees an intimate phase exchange. The construction of the separation chamber safely prevents carrying over of liquid particles into the vapour phase. Sample take-off from the vapour (6) and liquid (7) phases into vacuum and pressure-tight removable receivers is effected via manually controlled solenoid valves (5). Both samples may also be directly extracted from the circulation streams (8-9) by means of a gas-tight syringe. This feature is especially useful for the subsequent gaschromatographical evaluation.

A special advantage of the apparatus is that in case of mixtures with miscibility gaps, a direct gaseous sampling of vapour phase is possible by means of another sampling port (10).

The apparatus is designed for a charge of approx. 35 ml and can be operated under vacuum, at atmospheric or over-pressure and at temperatures of up to approx. 250°C. The electronic control is accommodated in a stylish housing.

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Scientists and chemical engineers in industry, research and teaching still measure phase-equilibria, despite the improvements in theoretical calculations.

For that reason, the FISCHER-phase equilibrium apparatus is part of the standard equipment in chemical engineering laboratories.

The vapour-liquid equilibria of binary and multi-component mixtures constitute an important basis for the determination of the theoretical plates necessary for separation. Although data is known from literature, there is an increasing tendency towards re-determination of equilibria data by experiment because the thermodynamic parameters such as pressure, temperature and concentrations are different or inadequate.

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Fig. 1
Relation between an equilibrium curve and a plate column with ideally working plates.
Graphical determination of the theoretical number of plates by drawing the steps in the McCabe Thiele diagram.